

Brilliant Violet 650 anti-human CD10

Mouse IgG1, kappa
HI10a



BioLegend
From Revvity

100 tests

Artikel-Nr.: B312242 | **BioLegend** | Hersteller-Nr.: 312242

364,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: HI10a
Format: BV650
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 Tests

Brilliant Violet 650 anti-human CD10

Clone: HI10a

Isotype: Mouse IgG1, kappa

Other_Names: Common acute lymphoblastic leukemia antigen (CALLA), Enkephalinase, Neutral endopeptidase, Neprilysin

Size: 100 tests

Reactivity: Human, Cynomolgus, Rhesus

Format: BV650

App_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=23965>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>