

Brilliant Violet 650 anti-human CD158e1 (KIR3D...

Mouse IgG1, kappa
DX9



BioLegend
From Revvity

100 tests

Artikel-Nr.: B312730 | **BioLegend** | Hersteller-Nr.: 312730

356,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: DX9
Format: BV650
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 Tests

Brilliant Violet 650 anti-human CD158e1 (KIR3DL1, NKB1)

Clone: DX9

Isotype: Mouse IgG1, kappa

Other_Names: NKAT3, NKB1, KIR-NKB1, p70

Size: 100 tests

Reactivity: Human

Format: BV650

App_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=23894>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>