

Brilliant Violet 421 anti-human CD268 (BAFF-R)

Mouse IgG1, kappa
11C1



BioLegend
From Revvity

100 tests

Artikel-Nr.: B316932 | BioLegend | Hersteller-Nr.: 316932

348,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: 11C1
Format: BV421
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 Tests

Brilliant Violet 421 anti-human CD268 (BAFF-R)

Clone: 11C1

Isotype: Mouse IgG1, kappa

Other_Names: TNFRSF13C, BAFF-R, BAFFR, BR3, BAFF Receptor

Size: 100 tests

Reactivity: Human

Format: BV421

App_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=23900>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>