

TotalSeq-B0593 anti-human CD203c (E-NPP3)

Mouse IgG1, kappa
NP4D6



BioLegend
From Revvity

10 µg
Artikel-Nr.: B324633 | **BioLegend** | Hersteller-Nr.: 324633

296,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: NP4D6
Format: TOTALSEQ_B
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 10 µg

TotalSeq-B0593 anti-human CD203c (E-NPP3)

Clone: NP4D6
Isotype: Mouse IgG1, kappa
Other_Names: E-NPP3, ENPP3
Size: 10 µg

Reactivity: Human
Format: TOTALSEQ_B
App_Abbrev: FC, ICFC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=23825>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>