

Brilliant Violet 421 anti-human CD85g (ILT7)

Mouse IgG1, kappa
17G10.2



BioLegend
From Revvity

100 tests

Artikel-Nr.: B326420 | BioLegend | Hersteller-Nr.: 326420

316,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: 17G10.2
Format: BV421
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 Tests

Brilliant Violet 421 anti-human CD85g (ILT7)

Clone: 17G10.2

Isotype: Mouse IgG1, kappa

Other_Names: CD85g, immunoglobulin-like transcript 7, LILRA4

Size: 100 tests

Reactivity: Human

Format: BV421

App_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=23940>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>