

TotalSeq-Bn1316 anti-human HLA-DR

Mouse IgG2b, kappa
LN3



BioLegend
From Revvity

10 µg
Artikel-Nr.: B327027 | **BioLegend** | Hersteller-Nr.: 327027

296,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: LN3
Format: TOTALSEQ_BN
Isotype: Mouse IgG2b, kappa
Species of Origin: Mouse
Verpackung: 10 µg

TotalSeq-Bn1316 anti-human HLA-DR

Clone: LN3

Isotype: Mouse IgG2b, kappa

Other_Names: Major Histocompatibility Class II, MHC class II

Size: 10 µg

Reactivity: Human

Format: TOTALSEQ_BN

App_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=23276>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>