

PE/Cyanine5 anti-human CD39

Mouse IgG1, kappa
A1



BioLegend
From Revvity

100 tests

Artikel-Nr.: B328248 | **BioLegend** | Hersteller-Nr.: 328248

272,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: A1
Format: PE_CYANINE5
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 Tests

PE/Cyanine5 anti-human CD39

Clone: A1
Isotype: Mouse IgG1, kappa
Other_Names: gp80, E-ATPDase, NTPDase-1, ecto-apyrase, Ec3.6.1.5
Size: 100 tests

Reactivity: Human, Cynomolgus, Rhesus
Format: PE_CYANINE5
App_Abbrev: FC, ICFC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=23834>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>