

Human CD3/CD28 T Cell Activation Beads



BioLegend
From Revvity

2.5 mL

Artikel-Nr.: B422604 | BioLegend | Hersteller-Nr.: 422604

492,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: Activ, Agonist, Cell Culture, FC, Stim
Format: SOLUTION
Species of Origin: N/A
Verpackung: 2.5 ml

Human CD3/CD28 T Cell Activation Beads

Clone:

Isotype:

Other_Names: T Cell Stimulation Beads, Human T Cell Stimulation Beads, Human CD3 CD28 Beads, CD3, CD28

Size: 2.5 mL

Reactivity: Human

Format: SOLUTION

App_Abbrev: FC, Activ, Stim, Cell Culture, Agonist

Species of Origin: N/A

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=23351>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>