

## Caretta 96-Well PCR Platten, farblos, weiß, barcod

PC Full Skirt, Low Profile, Barcode, Eckschnitt H1  
DNA-, DNase-, RNasefrei



Packung à 100 Stück  
Artikel-Nr.: 711812N | Biozym

**764,00 €\***

\*zzgl. MwSt. [zzgl. Versandkosten](#)

### Beschreibung

**Cyclertyp:** Analytic Jena / Biometra-FlexCycler2, Analytic Jena / Biometra-T1, Analytic Jena / Biometra-TAdvanced 96, Analytic Jena / Biometra-TAdvanced 96G, Analytic Jena / Biometra-TAdvanced 96S, Analytic Jena / Biometra-TAdvanced 96SG, Analytic Jena / Biometra-TGradient, Analytic Jena / Biometra-TOne, Analytic Jena / Biometra-TProfessional Serie, Analytic Jena / Biometra-qTOWERiris (touch), Bio-Rad / MJ Research-C1000 (Touch), Bio-Rad / MJ Research-CFX96 Touch Deep Well, Bio-Rad / MJ Research-DNA Engine Tetrad (2), Bio-Rad / MJ Research-DNA Engine, Bio-Rad / MJ Research-DNA Engine Dyad, Bio-Rad / MJ Research-DNA Engine Dyad Disciple, Bio-Rad / MJ Research-PC Tempo 96 / Deep Well, Bio-Rad / MJ Research-PTC 200, Bio-Rad / MJ Research-S1000, Bioer-GeneTouch 96, Bioer-GeneTouch Plus 96, Bioer-LifeTouch, Bioer-Lifeeco, Eppendorf-Mastercycler ep (gradient), Eppendorf-Mastercycler ep pro (S), Eppendorf-Mastercycler ep realplex, Eppendorf-Mastercycler nexus (gradient), Eppendorf-Mastercycler nexus GSX1, Eppendorf-Mastercycler nexus SX1, VWR-Collection UNO96

**Eckschnitt:** H1

**Eigenschaften:** mit Barcode, rigid frame

**Farbe:** Rahmen: farblos, weiß

**Format:** full skirted

**Kategorie:** 96-Well Platte

**Produktlinie:** Caretta

**Profil:** Low Profile

**Verpackung:** Packung à 100 Stück

**Volumen:** 0.1 ml