


96 Well PCR Plate, clear

High Profile, cut A12
Semi-Skirt



Packung à 50 Stück

Artikel-Nr.: AF4TI-0760 | Azenta | Hersteller-Nr.: 4TI-0760

183,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

Eckanschnitt: A12

Farbe: farblos

Format: semi skirted

Kategorie: 96-Well Platte

Produktlinie: Azenta

Profil: High Profile

Verpackung: Packung à 50 Stück

Volumen: 0.2 ml

Diese 96-Well PCR Platte ist kompatibel mit folgenden Thermocyclern, qPCR Cyclern oder Sequencern:

ABI®/LIFE TECHNOLOGIES/THERMO FISHER SCIENTIFIC

- Thermal Cyclers (96 well standard block): Veriti | Proflex | Simplicamp | GeneAmp® 2700/2720/9600/9700
- qPCR Cyclers (96 well standard block): 7000 | 7300 | 7500 | 7700 | 7900 HT | QuantStudio™ 3/5/6/7/12K | ViiA7™
- Sequencers (96 well standard block): 3100 | 3130XL | 3500 | 3500XL | 3730 | 3730XL

AGILENT/STRATAGENE

- qPCR Cyclers (96 well block): Mx3000P™ | Mx3005P™ | Mx4000™

ANALYTIK JENA/BIOMETRA

- Thermal Cyclers (96 well block): Flexcycler2 | T1 Thermocycler | Tgradient | Tone | Tadvanced | TProfessional (Standard/Basic) Gradient/XL | Trotbot 96 | SpeedCycler2 (SP, SPR)
- qPCR Cyclers (96 well block): qTOWER3/G/touch | Topical

BIOER TECHNOLOGIES

- Thermal Cyclers (96 well block): Gene Touch 96

BIO-RAD

- Thermal Cyclers (96 well block): C1000 Touch | S1000 | iCycler™ | MyCycler™ | T100
- qPCR Cyclers (96 well block): MyiQ™ | iCycler IQ™/IQ 4/IQ 5

BIO-RAD MJ RESEARCH

- Thermal Cyclers (96 well block): PCT100™/200™/220™/221™/225™/240™
- qPCR Cyclers (96 well block): Chromo4™

CORBETT RESEARCH

- Thermal Cyclers (96 well block): (Qiagen) Palm Cycler

EPPENDORF

- Thermal Cyclers (96 well block): MasterCycler® ep/ep gradient/Pro/Pro S/nexus/nexus gradient/nexus SX1/nexus GSX1
- qPCR Cyclers (96 well block): MasterCycler™ ep realplex

PEQLAB/VWR

- Thermal Cyclers (96 well block): peqSTAR 96X

SENSOQUEST

- Thermal Cyclers (96 well block): Labcycler

TAKARA

- Thermal Cyclers (96 well block): Dice touch | Gradient

TECHNE

- Thermal Cyclers (96 well block): Prime | PrimeG | Prime Elite | Prime Elite Satellite | PCRmax Alpha cyclers 1/2/4 | TC412 | TC512 | Genius | Genius Quad | Touchgene | Touchgene Gradient | Flexigene