

8 Well PCR Tube Strips, clear

High Profile



Packung à 125 Strips
Artikel-Nr.: AF4TI-0781 | Azenta | Hersteller-Nr.: 4TI-0781

80,50 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

Farbe: farblos
Kategorie: 8er-Streifen
Produktlinie: Azenta
Profil: High Profile
Verpackung: Packung à 125 Strips
Volumen: 0.2 ml

8 well PCR tube strips passend für:

ABI®/LIFE TECHNOLOGIES/THERMO FISHER SCIENTIFIC

- Thermal Cyclers (96 well standard block): Veriti | Proflex | Simpliamp | GeneAmp® 2700/2720/9600/9700
- qPCR Cyclers (96 well standard block): 7000 | 7300 | 7500 | 7700 | 7900 HT | QuantStudio™ 3/5/6/7/12K | ViiA7™

AGILENT/STRATAGENE

- Thermal Cyclers (96 well block): Surecycler 8800
- qPCR Cyclers (96 well block): Mx3000P™ | Mx3005P™ | Mx4000™

ANALYTIK JENA/BIOMETRA

- Thermal Cyclers (strips only): TRIO | Tpersonal | T3 Thermocycler
- Thermal Cyclers (96 well block): Flexcycler2 | T1 Thermocycler | Tgradient | Tone | Tadvanced | TProfessional (Standard/Basic) Gradient/XL | Trobot 96 | SpeedCycler2 (SP, SPR)
- qPCR Cyclers (96 well block): qTOWER3/G/touch | Toptical

BIOER TECHNOLOGIES

- Thermal Cyclers (strips only): GeneQ
- Thermal Cyclers (96 well block): Gene Touch 96

BIO-RAD

- Thermal Cyclers (strips only): Genecycler
- Thermal Cyclers (96 well block): C1000 Touch | S1000 | iCycler™ | MyCycler™ | T100

BIO-RAD MJ RESEARCH

- Thermal Cyclers (strips only): Mini Gradient
- Thermal Cyclers (96 well block): PCT100™/200™/220™/221™/225™/240™
- qPCR Cyclers (96 well block): Chromo4™

CORBETT RESEARCH

- Thermal Cyclers (96 well block): (Qiagen) Palm Cyclers

EPPENDORF

- Thermal Cyclers (96 well block): MasterCycler® ep/ep gradient/Pro/Pro S/nexus/nexus gradient/nexus SX1/nexus GSX1 | MasterCycler™ nexus X2/GX2/GX2e/ X2e
- qPCR Cyclers (96 well block): MasterCycler™ ep realplex

PEQLAB/VWR

- Thermal Cyclers (strips only): peqSTAR XS | peqSTAR 2X
- Thermal Cyclers (96 well block): peqSTAR 96X

SENSOQUEST

- Thermal Cyclers (96 well block): Labcycler

TAKARA

- Thermal Cyclers (96 well block): Dice touch | Gradient

TECHNE

- Thermal Cyclers (strips only): 3Prime | 3PrimeG | 3PrimeX
- Thermal Cyclers (96 well block): Prime | PrimeG | Prime Elite | Prime Elite Satellite | PCRmax Alpha cyclers 1/2/4 | TC412 | TC512 | Genius | Genius Quad | Touchgene | Touchgene Gradient | Flexigene