


FrameStar 96 PCR Plate, clear Wells, blue Frame

Low Profile, cut A12
Semi-Skirt, breakable vertically



Packung à 50 Stück
Artikel-Nr.: AF4TI-1200-B | Azenta | Hersteller-Nr.: 4TI-1200/B

302,50 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

Eckanschnitt: A12
Eigenschaften: breakable, rigid frame
Farbe: blau
Farbe Wells: farblos
Format: semi skirted
Kategorie: 96-Well Platte
Produktlinie: FrameStar
Profil: Low Profile
Verpackung: Packung à 50 Stück
Volumen: 0.1 ml

Diese 96-Well FrameStar Platte ist kompatibel mit folgenden Thermocyclern, qPCR Cyclern oder Sequencern

ABI®/LIFE TECHNOLOGIES/THERMO FISHER SCIENTIFIC

- GeneAmp® 9800 FAST | Verity FAST
- StepOne
- 7500 FAST | 7900 HT FAST | ViiA7™

AGILENT/STRATAGENE

- Robocycler Gradient
- AriaMx

BIO RAD

- C1000 Touch | S1000
- CFX96 Touch | CFX96 Touch Deep Well | CFX connect

BIO-RAD MJ RESEARCH

- PCT100™/200™/220™/221™/225™/ 240™
- MiniOpticon (strips only)
- Opticon2™

EPPENDORF

- MasterCycler® ep/ep gradient/Pro/Pro S/nexus/nexus gradient/nexus SX1/nexus GSX1 | MasterCycler™ nexus X2/GX2/GX2e/ X2e (strips only)
- MasterCycler™ ep realplex

PEQLAB/VWR

- peqSTAR XS | peqSTAR 2X (strips only)
- peqSTAR 96X

SENSOQUEST

- Labcycler

TAKARA

- Dice touch | Gradient

TECHNE

- TC412 | TC512 | Genius | Genius Quad | Touchgene | Touchgene Gradient | Flexigene
- Quantica