

## FrameStar 96 PCR Plate, clear Wells, red Frame

Low Profile, cut A12  
Semi-Skirt, breakable vertically



SCIENCE IS OUR BUSINESS



Packung à 50 Stück

Artikel-Nr.: AF4TI-1200-R | Azenta | Hersteller-Nr.: 4TI-1200/R

**302,50 €\***

\*zzgl. MwSt. zzgl. [Versandkosten](#)

### Beschreibung

**Eckanschnitt:** A12  
**Eigenschaften:** breakable, rigid frame  
**Farbe:** rot  
**Farbe Wells:** farblos  
**Format:** semi skirted  
**Kategorie:** 96-Well Platte  
**Produktlinie:** FrameStar  
**Profil:** Low Profile  
**Verpackung:** Packung à 50 Stück  
**Volumen:** 0.1 ml

**Diese 96-Well FrameStar Platte ist kompatibel mit folgenden Thermocyclern, qPCR Cyclern oder Sequencern**

#### ABI®/LIFE TECHNOLOGIES/THERMO FISHER SCIENTIFIC

- GeneAmp® 9800 FAST | Verity FAST
- StepOne
- 7500 FAST | 7900 HT FAST | ViiA7™

#### AGILENT/STRATAGENE

- Robocycler Gradient
- AriaMx

## BIO RAD

- C1000 Touch | S1000
- CFX96 Touch | CFX96 Touch Deep Well | CFX connect

## BIO-RAD MJ RESEARCH

- PCT100™/200™/220™/221™/225™/ 240™
- MiniOpticon (strips only)
- Opticon2™

## EPPENDORF

- MasterCycler® ep/ep gradient/Pro/Pro S/nexus/nexus gradient/nexus SX1/nexus GSX1 | MasterCycler™ nexus X2/GX2/GX2e/ X2e (strips only)
- MasterCycler™ ep realplex

## PEQLAB/VWR

- peqSTAR XS | peqSTAR 2X (strips only)
- peqSTAR 96X

## SENSOQUEST

- Labcycler

## TAKARA

- Dice touch | Gradient

## TECHNE

- TC412 | TC512 | Genius | Genius Quad | Touchgene | Touchgene Gradient | Flexigene
- Quantica