

384 Well PCR Plate, clear

cut A24
Skirt



Packung à 50 Stück

Artikel-Nr.: AF4TI-1384 | Azenta | Hersteller-Nr.: 4TI-1384

196,50 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

Eckanschnitt: A24

Farbe: farblos

Format: full skirted

Kategorie: 384-Well Platte

Produktlinie: Azenta

Verpackung: Packung à 50 Stück

Volumen: 30 µl

Diese 384-Well PCR Platte ist kompatibel mit folgenden Thermocyclern, qPCR Cyclern oder Sequencern

ABI®/LIFE TECHNOLOGIES/THERMO FISHER SCIENTIFIC

- Thermal Cyclers (384 well block): GeneAmp® 9700 | Verity | Proflex | Multiblock system
- qPCR Cyclers (384 well block): QuantStudio™ 5/6/7/12K | ViiA7™ | 7900 HT FAST
- Sequencers (384 well block): 3100 | 3130XL | 3500 | 3500XL | 3730 | 3730XL

AGILENT/STRATAGENE

- Thermal Cyclers (384 well block): SureCycler 8800

ANALYTIK JENA/BIOMETRA

- Thermal Cyclers (384 well block): Flexcycler2 | T1 Thermocycler | Tadvanced | TProfessional | Trobot 96
- qPCR Cyclers (384 well block): qTOWER3 84 / 84G

BIOER TECHNOLOGIES

- Thermal Cyclers (384 well block): Gene Touch 384

BIO-RAD

- Thermal Cyclers (384 well block): C1000 Touch | S1000
- qPCR Cyclers (384 well block): CFX384 Touch

BIO-RAD MJ RESEARCH

- Thermal Cyclers (384 well block): PCT200™/220™/221™/225™/240™

CORBETT RESEARCH

- Thermal Cyclers (384 well block): (Qiagen) Palm Cyclyer 384

EPPENDORF

- Thermal Cyclers (384 well block): MasterCycler® ep 384/Pro 384

GE HEALTHCARE/AMERSHAM

- Sequencers (384 well block): MegaBACE™ 4000

PEQLAB/VWR

- Thermal Cyclers (384 well block): peqSTAR 384X

SENSOQUEST

- Thermal Cyclers (384 well block): Labcycler

TECHNE

- Thermal Cyclers (384 well block): Prime | PrimeG | Prime Elite | Prime Elite Satellite | PCRmax Alpha cycler 1/2/4 | TC412 | TC512 | Genius | Genius Quad | Flexigene