

IntelliXcap Automated Septum Cap Decapper ReCapper

96-format, for validated sample storage tubes and TPE caps, purge module, vibration chute



1 Stück

Artikel-Nr.: AX46-5030 | Azenta | Hersteller-Nr.: 46-5030

54.480,50 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

Verpackung: 1 Stück

IntelliXcap Septum Caps, automated Septum Cap Decapper Re-Capper, 96-format, compatible with validated sample storage tubes and TPE caps, fitted with a purge module that allows an inert gas of choice to be layered over the samples and vibration chute for cap disposal