

IntelliXcap Septum Caps Adaptor Kit for Azenta

for 0.9 ml Int. Thread (66-62321-L), 0.48 ml Int.
0.3ml Int. Thread, TPE Septum Cap, Int. Thread



1 Stück

Artikel-Nr.: AX46-5050 | Azenta | Hersteller-Nr.: 46-5050

1.534,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

Verpackung: 1 Stück

IntelliXcap Septum Caps Adaptor Kit, for Azenta, compatible with 0.9ml Internal Thread (66-62321-L), 0.48ml Internal Thread Jacket (67-0753-01), 0.3ml Internal Thread Jacket (66-62325), TPE Septum Cap, Internal Thread (65-73000)