

## Rack, 96-format, HighBase, standard profile

non-locking lid, for 1.0 ml Tubes with Septum TPE  
Caps/Screw Caps, 1.4 ml Tubes with Septum Caps



SCIENCE IS OUR BUSINESS



Packung à 10 Racks

Artikel-Nr.: AX66-51001 | Azenta | Hersteller-Nr.: 66-51001

**180,00 €\***

\*zzgl. MwSt. zzgl. Versandkosten

### Beschreibung

<b>Brand:</b>	FluidX
<b>Codierung:</b>	ohne Codierung
<b>Format:</b>	HighBase Rack
<b>Passend für Rackformat:</b>	SBS/96
<b>Produktart:</b>	Rack leer
<b>Verpackung:</b>	Packung à 10 Racks
<b>Volumen:</b>	1.0 ml, 1.4 ml