

Purified anti-human CD90 (Thy1)

Mouse IgG1, kappa
S22015B



BioLegend
From Revvity

100 µg

Artikel-Nr.: B389802 | BioLegend | Hersteller-Nr.: 389802

136,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC
Clone: S22015B
Format: PURE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 µg

Purified anti-human CD90 (Thy1)

Clone: S22015B
Isotype: Mouse IgG1, kappa
Other_Names: Thy-1, Thy1
Size: 100 µg

Reactivity: Human
Format: PURE
App_Abbrev: FC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=24607>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>