

Purified anti-human CD138 (Syndecan-1)

Rat IgG2b, kappa
S22003F



BioLegend
From Revvity

100 µg

Artikel-Nr.: B390206 | BioLegend | Hersteller-Nr.: 390206

108,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC
Clone: S22003F
Format: PURE
Isotype: Rat IgG2b, kappa
Species of Origin: Rat
Verpackung: 100 µg

Purified anti-human CD138 (Syndecan-1)

Clone: S22003F
Isotype: Rat IgG2b, kappa
Other_Names: Syndecan-1, B-B4
Size: 100 µg

Reactivity: Human
Format: PURE
App_Abbrev: FC
Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=24759>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>