

2X RT-RPA exo Master Mix, lyophilized



8 rxn of 50 µl

Artikel-Nr.: 351600400 | biotechrabbit | Hersteller-Nr.: BC1600400

44,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

Format: Lyophilisiert

Kategorie: Recombinase Polymerase Reaktion

Produkttyp: Kit

Verpackung: 8 rxn of 50 µl

Storage Temperature: 4°C or RT

Shipping Temperature: Ambient, Coolpack or Dry Ice

Minimal remaining shelf life (months): 6