

LYO-ready T4 Gene 32 Protein, 10 mg/ml



500 µg

Artikel-Nr.: 351500501 | biotechrabbit | Hersteller-Nr.: BR1500501

288,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

Format: Lyo-ready
Kategorie: Recombinase Polymerase Reaktion
Produkttyp: Komponente
Verpackung: 500 µg

Storage Temperature: -20°C
Shipping Temperature: Coolpack or Dry Ice
Minimal remaining shelf life (months): 6