

## LYO-ready T4 Gene 32 Protein, 10 mg/ml



500 µg

Artikel-Nr.: 351500501 | biotechrabbit | Hersteller-Nr.: BR1500501

**288,00 €\***

\*zzgl. MwSt. zzgl. Versandkosten

### Beschreibung

**Format:** Lyo-ready  
**Kategorie:** Recombinase Polymerase Reaktion  
**Produkttyp:** Komponente  
**Verpackung:** 500 µg

Storage Temperature: -20°C

Shipping Temperature: Coolpack or Dry Ice

Minimal remaining shelf life (months): 6