

## Spark PLUS UV395 anti-human CD326 (EpCAM)

Mouse IgG2b, kappa  
9C4



BioLegend  
From Revvity

25 tests

Artikel-Nr.: B324259 | BioLegend | Hersteller-Nr.: 324259

**177,00 €\***

\*zzgl. MwSt. zzgl. Versandkosten

### Beschreibung

**App\_Abbrev:** FC  
**Clone:** 9C4  
**Format:** SPKPLUSUV\_395  
**Isotype:** Mouse IgG2b, kappa  
**Species of Origin:** Mouse  
**Verpackung:** 25 tests

Spark PLUS UV395 anti-human CD326 (EpCAM)

Clone: 9C4

Isotype: Mouse IgG2b, kappa

Other\_Names: Ep-CAM, tumor associated calcium signal transducer 1, epithelial cell surface antigen, epithelial glycoprotein 2, EGP2, adenocarcinoma associated antigen, TROP1

Size: 25 tests

Reactivity: Human

Format: SPKPLUSUV\_395

App\_Abbrev: FC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=26729>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>